# FREEPRINT® splint

WAVELENGTH

**385 nm !**

fast building speed

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Resin for dental 3D printing

## INDICATION

Drilling and X-ray templates for implant dentistry, occlusal splints, fixation and transfer keys

## PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

## INTENDED USERS

Dentist, dental technician

## SUITABLE FOR THE FOLLOWING DLP-PRINTER

Asiga 385 nm		MiiCraft 385 nm		Rapidshape 385 nm		W2P 385 nm
385 nm	405 nm			385 nm	405 nm	
MAX UV	MAX 405	MiiCraft y-Serie		DII-Serie	D-Serie	Solflex-Serie
Pro2 UV	Pro2 405					
Pico2 UV	Pico2 405					

Printers may only be operated using material parameters authorised by DETAX!

## PROCESSING

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe processdescription on page 2).
- ▶ After storage, the material in the bottle should be shaken intensively and homogenized with a bottle roller before use.
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ **Maximum curing depth\* at direct post-exposure:** clear: 3,5 mm  
\*In case of large objects and exposure on both sides, the material thickness can be up to 7 mm (Example FREEPRINT® splint – with a curing depth of 3,5 mm).  
Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ▶ Seal surface with FREEFORM® coat.
- ▶ Processing temperature 23 °C ± 2 °C.

## SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritation to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

## CLEANING

### By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

### By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards. Cleaning tablets can be used following the manufacturer's instructions.

### Note

- ▶ Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- ▶ Do not use any acids or solvents.

## NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.
- ▶ Please follow the instructions on the safety data sheet!

### For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.de and to the competent authority of the Member State in which the user and/or patient is established.

## STORAGE

FREEPRINT® splint is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

## CONTRAINDICATION

Contains (meth)acrylics and phosphine oxide.

Some ingredients of FREEPRINT® splint may cause allergic reactions in predisposed persons. In such cases refrain from using the product. FREEPRINT® splint only insert intraorally in completely polymerised state.

## ADVERSE EFFECTS

Product may cause allergic reactions.

## DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

### Symbols:

 Medical Device

### Storage:

### Application:

At 23 °C ± 2 °C

# Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol purity  $\geq 98\%$  for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

## Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C to remove the solvent residues from the cleaning process.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit (e.g. Otoflash G171) with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

## Surface processing

Varnish the surface with FREEFORM® coat

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GmbH & Co. KG